QUICK START GUIDE

C32.60.PM.05

Radwag

Distributed by AKRIBIS

IMMB-09-02-05-22-EN-PL

Start up Guide

Szybki start


Precision Balances

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- PM C32

To view a complete User Manual go to the website or scan QR code:

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






PM C32 (max: 10-50 kg)

PM C32 (max ≥ 60 kg)


Keep the packaging should a warranty claim or service be required.

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3. RECOMMENDATIONS AND WARNINGS

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
Operate the device in a room where the temperature ranges between 10–40 °C (50–104 °F) and where the relative humidity is below 80%. Temperature must change slowly and gradually.


Using the anti-vibration table is recommended in order to enable the most stable and repeatable measurement.

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
Place the balance away from heat sources. Avoid exposing the balance to the sunlight.

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
Avoid exposing the balance to a magnetic field. Do not weigh magnetic substances.

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
Avoid air drafts and air movements at the workstation.

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
Make sure that the balance is placed on an even surface.

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Do not place the balance on an unstable surface exposed to shocks and vibrations.

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Not intended for usage within EX zones. Not intended for weighing explosive or flammable substances.

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
4. ACTIVITIES TO BE DONE PRIOR TO OPERATION

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


Check whether the grounding spring takes its intended location. Make sure that the grounding spring juts slightly out of the hole.

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

PM C32 (max ≥ 60 kg)

6. PREPARING FOR WORK

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
Place the device at the workstation, wait until it reaches room temperature.

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6.1. Temperature stabilization

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
1

Prior to switching the device on, let it reach room temperature (estimated stabilisation time: 1-8 hours).

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
2

Connect the power adapter to DC connector.

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3


Connect the power adapter to the mains.

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6.2. Balance levelling

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

PM C32 (max: 10 - 50 kg)

To level the balance, turn its feet and keep turning the feet until an air bubble takes the central position.

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
6.3. Switching the balance on

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To switch the device on, press ON / OFF key.

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





8.2. Zeroing

Zerowanie

Make sure that the weighing pan is empty and press **+0+** key.

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

8.3. Taring

Tarowanie

With loaded weighing pan: upon stabilization of the weighing result, press **T+** key.

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


Using keys **↑** **↓** select working mode and press **✓** key to confirm.

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8.5. Units selection

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Press the **UNIT** key to change the weight unit.

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The balance enables selecting the following weighing units:

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9. SETTINGS

USTAWIENIA

Some balance settings are accessible for Administrator exclusively.


Prior balance parameters setup, log in as the Administrator.

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


Users list is displayed. Using  keys, select "Admin" option and press  key to confirm.

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Using numeric keyboard enter "1111" password and press  key to confirm.

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You have logged as the administrator.

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9.2. User settings

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
Prior balance users setup, log in as Administrator (point 9.1.).

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
1

Press the key under the  icon to enter databases.

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
2

Using  keys, go to "Users" and press  key to confirm.

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


5

Using keyboard enter user name and press  key to confirm.

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Password: Enter user password

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Access levels: Set respective access level (user, advanced, admin)

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New	
Name	John Smith
Code	
Password	*****
Permissions	User

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
Prior balance proximity sensors setup, log in as Administrator (point 9.1.).

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
You can assign the proximity sensors with different functions for each working mode (the sensors will trigger different operation for each working mode). The diagram presents sensors setup for the "Weighing" mode.

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

Using  keys, go to "Working Modes" and press  key to confirm.

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



11. DEVICE CLEANING

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
11.1. General recommendations

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
Prior to cleaning it is necessary to disconnect the device from the mains! Remove the plug from the outlet, disconnect all cables.

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Use soft cloths made of microfiber, natural fiber or man-made fiber. Avoid using abrasive cloths or cloths that might scratch the surface.

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Apply the cleanser onto the cloth first. Avoid applying the cleanser directly onto the device.

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Do not use aggressive cleaning agents (e.g. solvents, chlorine preparations, corrosive chemicals, bleach).

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Neither use substances of sharp structure nor intended for scrubbing.

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Make sure that neither dust nor liquid gets inside the weighing system (device inside).

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