QUICK START GUIDE

PUE C32

Radwag

Distributed by AKRIBIS

ITKB-09-02-05-22-EN-PL

Start up Guide

Szybki start

Indicator

Miernik

- PUE C32

To view the full user manual, go to the website or scan the QR code:

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Indicator
Miernik
x1


Power supply*
Zasilacz sieciowy*
x1

**The type of plug may vary depending on the country.*

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
2. NOTES AND WARNINGS

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
Operate the indicator in a room free from dust where the temperature ranges -10 – 40 °C ($+14$ – 104 °F) and where the relative humidity is below 80%. Potential temperature change shall occur gradually and slowly. Place the scale on a table or a wall bracket ensuring stability. For the most stable and repeatable results it is recommended to use an anti-vibration table.

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
Place the scale away from heat sources. Avoid exposing the scale to the sunlight.

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Avoid exposing the scale to the influence of a magnetic field. Do not weigh magnetic substances.

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Not intended for use within EX zones. Not intended for weighing explosive or flammable substances.

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3. PREPARING FOR WORK

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
Place the indicator on the workstation and leave it until it reaches ambient temperature.

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

Prior to switching the indicator on let it reach room temperature (estimated stabilization time: 8 hours).

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Connect the power supply to a DC connector located on the side of the indicator.

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
Connect the power supply to the mains.

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3.2. Indicator Start-Up

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Press ON/OFF key to start the indicator.

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The indicator is ready for operation!

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



5. SETTINGS

USTAWIENIA

Some scale settings are only available to the Administrator.


Prior to setting scale parameters log in as the Administrator.

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


A list of operators is displayed. Use  keys to select "Administrator" option and press  key to confirm.

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When the field for filling the password is opened, confirm the empty field by pressing  key.

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
Prior to the proximity sensors configuration, log in as the Administrator (see step 5.1.).

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
You can assign a function of your choice to each proximity sensor. In addition, you can assign a different function to the same sensor in different working modes. See the exemplary procedure for sensors configuration for the “Weighing” mode.

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
Press  key.

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
Using  keys select “Working modes” and press  key to confirm.

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Using keys  select „Weighing” mode and press  key to confirm.

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
Using  keys select ‘Proximity sensors’ and press  key to confirm.

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
5

Using  keys select right or left sensor and press  to confirm.

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6

Using  keys select a function to be assigned to the respective sensor. Press  key to confirm.

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
Follow the above procedure to configure proximity sensors for other working modes.

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
7. CLEANING AND MAINTENANCE

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
7.1. Notes and Warnings

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
Prior to cleaning it is necessary to disconnect the scale from the mains. Unplug the connector from the electrical socket.

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For cleaning use only soft cloths made of microfiber. Do not use any abrasive cloths or any material that may scratch the surface!

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Apply the cleanser onto the cloth first. Avoid applying the cleanser directly onto the device!

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Avoid using aggressive cleaning agents (e.g. solvents, cleansers containing chlorine, corrosive substances and bleach).

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Avoid using cleansers containing abrasive substances or scouring preparations.

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Make sure that no dust or liquid enters the sockets and the inside of the indicator.

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- Plastic elements can be cleaned with a soft cloth slightly moistened with a solution of water and a mild detergent (e.g. soap or dishwashing liquid).
- Dry cleaned elements using soft and dry cloth or use dust-free paper towel to absorb the remaining moisture.
- Do not install the components until they are completely dry.

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