

DATASHEET

DJ-04

Radwag

Distributed by AKRIBIS





More information on the website
radwag.com/en/info,w1,ODT

DJ-04 Antistatic Ionizer

WX-001-0097



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet


Physical parameters	
Neutralization method	Corona discharge, alternating emission of positive and negative ions
Distance form a sample	5 – 200 mm (recommended 5 – 50 mm)
Neutralization time	<5 s (typically <2 s)
Frame dimensions	268×110×295 mm
Packaging dimensions W x D x H	375×375×150 mm
Net weight	1.3 kg
Gross weight	2.5 kg
Electrical parameters	
Power supply	power supply: 100 – 240V AC 50/60Hz 1A Max; 15V DC 2.4A ionizer: 12-15V DC 2W
High voltage	5 kV DC, pulsed, peak value
Maximum short-circuit current	<2mA (output HV)
Environmental conditions	
Operating temperature	0 – +50 °C
Relative humidity	5% – 80% no condensation

Compatible with (Additional Fee)


XA 5Y.M.A.P Microbalance for pipettes calibration
AS X7 Analytical Balance
XA 5Y.M.A Microbalance
XA 5Y.A Analytical Balance
AS 5Y Analytical Balance
XA 5Y.F Analytical Balance for Filters

XA 5Y Analytical Balance
AS R2 PLUS Analytical Balance
XA 5Y.M Microbalance
AS R1 PLUS Analytical Balance
XA 5Y.M.A.S Microbalance for stents weighing


Device dimensions W x D x H



DJ-04



DJ-05



power supply